

# Tobacco Produces Two Nitrosamines: An Opinion

Alex John\*

Department of Pharmaceutical Sciences, Government Medical College, Kerala, India

## Introduction

Tobacco-specific nitrosamines (TSNA) are carcinogens that have been linked to tobacco, tobacco smoke, and nicotine-containing goods. N-nitrosornicotine (NNN) and 4-(methylnitrosamino)-1-(3-pyridyl)-1-butanone; 4-(N-methyl-N-nitrosamino)-1-(3-pyridyl)-1-butanone; nicotine-derived nitrosamine ketone (NNK) were both shown to cause lung cancer in mice in 1964 and 1980, respectively. TSNA can be found in trace amounts in freshly harvested tobacco, but their concentration varies depending on the type of tobacco and fertilisers used throughout the growing process. Nicotine nitrosation is the primary source of NNN and NNK, while nornicotine nitrosation can also produce NNN. This happens most often when tobacco and tobacco products are processed, cured, and stored. NNK and NNN are partially derived from the distillation of these nitrosamines, which are pre-formed in the tobacco, and NNK is also a result of thermal release of the matrix-bound form, while another fraction is pyrosynthesized by nitrosation of the respective alkaloid precursors, possibly with nitrogen oxides derived from the nitrate, which is present in high concentrations in some tobacco types. Sidestream smoke contains NNK and NNN, and their yields are similar to or two to five times higher than those found in mainstream smoking.

## Description

The recovery, sensitivity, and precision of the approach were all excellent. For NNN and NNK, the method's detection limits were 0.0108 ng/m<sup>3</sup> and 0.0136 ng/m<sup>3</sup>, respectively. The instrument's calibration range was 0.2–60 ng/mL. The method's determined lower working range limit (LWRL) for NNN was 0.126 ng/m<sup>3</sup>, while NNK's LWRL was 0.195 ng/m<sup>3</sup>. The technique was used to assess surrogate ambient aerosols produced by smoking machines. Because the retention of NNN and NNK in the bodies of consumers is not taken into consideration, this model produces a considerable overestimation of the likely influence of THS 2.2 and e-cigarettes on indoor air. As a result, the given values do not reflect a real-life situation. THS 2.2 had 0.0830 0.0153 ng/m<sup>3</sup> of NNN and 0.0653 0.0138 ng/m<sup>3</sup> of NNK in the surrogate environmental aerosols, 0.0561 0.0296 ng/m<sup>3</sup> of NNN in e-cigarettes, and 0.816 0.109 ng/m<sup>3</sup> of NNN and 4.13 1.04 ng/m<sup>3</sup> NNK in cigarettes. These values correspond to 10% of the measured ETS concentration for NNN in THS 2.2 environmental aerosols and 7% for e-cigarette aerosols. For NNK, the THS 2.2 value for the environmental aerosol was 2% of the ETS value.

The researchers devised and verified a sensitive and accurate method for analysing two TSNA (NNN and NNK) in indoor air. A unique strategy for the collection employing silica sorbent cartridges with simplified posterior sample preparation was used to reduce the steps during sample preparation as well as the efficient concentration of the target compounds before analysis. After that,

\*Address for Correspondence: Alex John, Department of Pharmaceutical Sciences, Government Medical College, Kerala, India; E-mail: john.1978@gmail.com

Copyright: © 2022 John A. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Received: 02 April, 2022, Manuscript No. JBABM-22-62181; Editor Assigned: 04 April, 2022, PreQC No. P- 62181; Reviewed: 16 April, 2022, QC No. Q- 62181; Revised: 21 April, 2022, Manuscript No. R- 62181; Published: 28 April, 2022, DOI: 10.37421/1948-593X.2022.14.323

the extracts were examined by isotope dilution LC-MS/MS. When compared to the traditional methods of trapping on Cambridge filters, our procedure enhanced sample preparation recoveries. Furthermore, compared to previous methodologies on TSNA studies in air, the LOD and LLOQ were reduced by one to two orders of magnitude [1-4].

Tobacco-specific nitrosamines (TSNA) are carcinogens that have been linked to tobacco, tobacco smoke, and nicotine-containing goods. N-nitrosornicotine (NNN) and 4-(methylnitrosamino)-1-(3-pyridyl)-1-butanone; 4-(N-methyl-N-nitrosamino)-1-(3-pyridyl)-1-butanone; nicotine-derived nitrosamine ketone (NNK) were both shown to cause lung cancer in mice in 1964 and 1980, respectively. TSNA can be found in trace amounts in freshly harvested tobacco, but their concentration varies depending on the type of tobacco and fertilisers used throughout the growing process. Nicotine nitrosation is the primary source of NNN and NNK, while nornicotine nitrosation can also produce NNN. This happens most often when tobacco and tobacco products are processed, cured, and stored. NNK and NNN are partially derived from the distillation of these nitrosamines, which are pre-formed in the tobacco, and NNK is also a result of thermal release of the matrix-bound form, while another fraction is pyrosynthesized by nitrosation of the respective alkaloid precursors, possibly with nitrogen oxides derived from the nitrate, which is present in high concentrations in some tobacco types. Sidestream smoke contains NNK and NNN, and their yields are similar to or two to five times higher than those found in mainstream smoking [5].

## Conclusion

The accuracy profile process allowed the method's performance to be evaluated as a function of the matrices. In the matrices under examination, the method's working ranges allowed quantification of the target analytes with a precision of 25%. The method's fitness for purpose for prospective comparative analyses of the ambient aerosol of THS 2.2 and e-cigarettes, as well as the ETS of cigarettes, was established by the validation results.

## Acknowledgement

None.

## Conflict of Interest

The author has no conflict of interest towards the article.

## References

1. Sleiman, Mohamad, Randy L. Maddalena, Lara A. Gundel and Hugo Destailats. "Rapid and sensitive gas chromatography-ion-trap tandem mass spectrometry method for the determination of tobacco-specific N-nitrosamines in secondhand smoke." *J Chromat* 1216 (2009): 7899-7905.
2. Lueso, María Gómez, Maya I. Mitova, Nicolas Mottier and Mathieu Schaller, et al. "Development and validation of a method for quantification of two tobacco-specific nitrosamines in indoor air." *J Chromat* 1580 (2018): 90-99.
3. Lee, Yoon-Seo, Ki-Hyun Kim, Sang Soo Lee and Richard JC Brown. "Analytical method for measurement of tobacco-specific nitrosamines in E-cigarette liquid and aerosol." *Appl Sci* 8 (2018): 2699.
4. Kim, Hyun-Ji and Ho-Sang Shin. "Determination of tobacco-specific nitrosamines in replacement liquids of electronic cigarettes by liquid chromatography-tandem mass spectrometry." *J Chromat* 1291 (2013): 48-55.

5. Forster, Mark, John McAughey, Krishna Prasad and Christopher Proctor. "Assessment of tobacco heating product THP1. 0. Part 4: Characterisation of indoor air quality and odour." *Regul Toxicol Pharmacol* 93 (2018): 34-51.

**How to cite this article:** John, Alex. "Tobacco Produces Two Nitrosamines: An Opinion." *J Bioanal Biomed* 14 (2022): 323.